

ECG Research Tool – User Manual

This tool is an EDUCATIONAL & RESEARCH-ONLY ECG interpretation system. It is NOT intended for clinical or diagnostic use. Predictions may be incorrect. For any health concerns, ALWAYS consult a qualified medical professional.

How to Use the Website

1. Open the Website

Run the server and open the link:
<http://127.0.0.1:5000/>

2. Choose Your Analysis Type

- PTB-XL (10-second ECG segment classification)
- MIT-BIH (single beat arrhythmia classification)

3. Upload an ECG CSV File

On the page, click “Choose File” and select your CSV file.
A sample CSV is available for testing.

4. Click “Analyze ECG”

The system will process your ECG, draw the waveform, and show predicted class probabilities.

5. View Results

You will see:

- ECG waveform graph
- Predicted class
- Confidence probabilities
- Safety disclaimer

Supported File Format

CSV file containing ECG values, one column per lead or one column for single-lead.

Disclaimer

This is a research project only and must NOT be used for any medical decision-making.